

## Ceti Magnum Infinity Binocular Compound Microscope

**Product Code:** 2629.0000M

**Brand:** Ceti

Ceti's Magnum Infinity Binocular Compound Microscope is used for Brightfield viewing in the field or in regions where power supplies are interrupted. Providing excellent quality images via the widefield FN20mm eyepiece, this microscope comes with four Infinity Corrected Plan Achromatic objectives.

It features a robust, T-shaped base, designed to offer maximum stability. Modern, ergonomic features such as a 30° inclined eyepiece tube, ensure a comfortable viewing position over a prolonged period.

This model also comes with a bright, white LED light source with intensity control.

### Features

- Sturdy base for comfortable viewing over long periods
- Wide-field 10X, FN20 eyepieces with dioptric adjustment
- Quintuple, reversed nosepiece with four Infinity Corrected Plan Achromatic objectives
- Wide range of accessories
- Bright, white LED illumination with intensity control
- Full Köhler illumination



# Specifications

<b>Optical Head</b>	Binocular, inclined at 30°
	Interpupillary Distance Adjustment (Siedentopf Type): 55 to 75mm
	Dioptric adjustment on one eyepiece
<b>Eyepieces</b>	10x/20mm wide-field
<b>Nosepiece</b>	Quintuple reversed and revolving with click stop
<b>Objectives</b>	4x, 10x, 40x, 100x Infinity Corrected Plan Achromatic
<b>Condenser</b>	Pre-condenser with field diaphragm
	Iris diaphragm
<b>Stage</b>	210 x 140mm with built-in mechanical stage and slide holder
	90 x 60mm movement, vernier scale
<b>Focusing Knobs</b>	Coaxial coarse and fine with stop
	Tension adjustment ring
<b>Mains Power</b>	220 - 240V/50 - 60Hz
<b>Illumination</b>	3W LED with intensity control
<b>Supplied With</b>	Blue and green filter, 2 spare fuses, immersion oil (10ml), soft rubber eyecups, dust cover, and instruction manual

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)