



LSFSF107UK 107L Laboratory Essential Freezer – Solid

LSFSF107UK

- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	495 x 620 x 830
Weight	33kg
Gross Capacity (Litres)	107 Litres
Net Capacity (Litres)	103 Litres
Unit External Width (mm)	495mm
Unit External Depth (mm)	620mm
Unit External Height (mm)	830mm
Product Weight (kg)	33kg
Packed Width (mm)	560mm
Packed Depth (mm)	850mm
Packed Height (mm)	655mm
Packed Weight (kg)	36.5kg
Country Of Origin	UK
Internal Dimensions (mm) (WxDxH)	430 x 700 x 440
Temperature Range	-18°C to -25°C
Colour	White/Green

External Finish	Painted Steel
Internal Finish	HIPS*
Type Of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged Reversible
Door Liner	Flat
Door Open Depth (mm)	1050mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes (2 Keys)
Adjustable Feet	1 Adjustable Front Foot
Porthole	No
Temperature Display On Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	No
Audible Alarms	Yes
Visual Alarms	Yes
Number Of Shelves	2
Shelf Size (mm) (WxDxH)	427 x 374 x 6
Shelf Material	Polythene Coated Wire
Number Of Baskets	1
Basket Material	Polythene Coated Wire
Basket Size (mm) (WxDxH)	395 x 180 x 150
Interior Light	No
Fan	No
Spark Free	Yes
Defrost Setting	Manual
Controller Type	Digital External
Number Of Probes	1

Probe Type	Air
Ambient Temperature Range	10°C to 32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply (V)	230V
Insulation Thickness (mm)	35mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System	No
Optional Digital Lock	No
Other Door Option	No



For more information contact us on
 email: enquiries@medlinescientific.com
 Tel: +44 (0)1865 400321