



LSFSF39DC-UK 39L Laboratory Freezer, Solid door

LSFSF39DC-UK

SPECIFICATION

Dimensions	475mm x 445mm x 565mm
Weight	18kg
Gross Capacity (Litres)	39L
Net Capacity (Litres)	35L
Unit External Width (mm)	475mm
Unit External Depth (mm)	445mm
Unit External Height (mm)	565mm
Product Weight (kg)	18kg
Packed Width (mm)	490mm
Packed Depth (mm)	530mm
Packed Height (mm)	620mm
Packed Weight (kg)	22.5kg
Country Of Origin	UK
Internal Dimensions (mm) (WxDxH)	370mm x 345mm x 404mm
Temperature Range	-18°C to -25°C
Colour	White/Green
External Finish	Painted Steel

Internal Finish	HIPS*
Type Of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	905mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet	2, Front
Porthole	Yes
Porthole Size	16mm
Porthole Location	Rear Left-hand Corner
Temperature Display On Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number Of Shelves	1
Number Of Positions	1
Shelf Size (mm) (WxD)	365mm x 190mm
Shelf Material	Polythene Coated Wire
Number Of Baskets	1
Basket Material	Polythene Coated Wire
Basket Size (mm) (WxDxH)	363mm x 150mm x 150mm
Interior Light	No
Fan	No
Spark Free	Yes
Defrost Setting	Manual
Controller Type	Digital/ External

Number Of Probes	1
Probe Type	Air
Ambient Temperature Range	+10°C to +32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply (V)	220-240V
Insulation Thickness (mm)	50mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System Available	No
Optional Digital Lock	No
Other Door Option	No



For more information contact us on
 email: enquiries@medlinescientific.com
 Tel: +44 (0)1865 400321