



LSFSF98BT-UK 98L Laboratory Plus Freezer, Solid Door

LSFSF98DC-UK

SPECIFICATION

Dimensions	595mm x 660mm x 850mm
Weight	37.25kg
Gross Capacity (Litres)	98 Litres
Net Capacity (Litres)	84 Litres
Unit External Width (mm)	595mm
Unit External Depth (mm)	660mm
Unit External Height (mm)	850mm
Product Weight (kg)	37.25kg
Packed Width (mm)	640mm
Packed Depth (mm)	660mm
Packed Height (mm)	910mm
Packed Weight (kg)	41.75kg
Country Of Origin	UK
Internal Dimensions (mm) (WxDxH)	420mm x 430mm x 640mm
Temperature Range	-18°C to -25°C
Colour	White/Green

External Finish	Painted Steel
Internal Finish	HIPS*
Type Of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1130mm
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet	4, All Adjustable
Porthole	Yes
Porthole Size	16mm
Porthole Location	Rear Bottom
Temperature Display On Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number Of Shelves	3
Number Of Shelf Positions	3
Shelf Size (mm) (WxDxH)	495mm x 325mm x 6mm
Shelf Material	Glass/HIPS
Drawer Size (mm) (WxDxH)	2 - 385mm x 395mm x 195mm, 1 - 385mm x 242mm x 205mm
Drawer Material	HIPS
Number of Baskets	0
Interior Light	No
Fan	No
Spark Free	Yes

Defrost Setting	Auto (Off Cycle)
Controller Type	Digital, External
Bluetooth	Yes
Probe Type	Air
Ambient Temperature Range	10°C to 32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply (V)	220-240V
Insulation Thickness (mm)	80mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System Available	No
Optional Digital Lock	No
Other Door Option	No



For more information contact us on
 email: enquiries@medlinescientific.com
 Tel: +44 (0)1865 400321