



LSFSR158BT-UK 158L Laboratory Plus Refrigerator – Solid Door

LSFSR158BT-UK

- ✓ Local monitoring and alarm history via Bluetooth
- ✓ 2°C to 10°C temperature range
- ✓ Spark free interior
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	595mm x 660mm x 850mm
Weight	33.25kg
Gross Capacity (Litres)	158 Litres
Net Capacity (Litres)	154 Litres
Unit External Width (mm)	595mm
Unit External Depth (mm)	660mm
Unit External Height (mm)	850mm
Product Weight (kg)	33.25kg
Packed Width (mm)	640mm
Packed Depth (mm)	660mm
Packed Height (mm)	910mm
Packed Weight (kg)	37.82kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	500mm x 450mm x 725mm

Temperature Range	+2°C to +10°C
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type Of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1150mm
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet	4, All Adjustable
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear Bottom
Temperature Display On Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number Of Shelves	4
Number Of Shelf Positions	5
Shelf Size (mm) (WxDxH)	495mm x 325mm x 6mm
Shelf Material	Polythene Coated Wire
Number Of Baskets	1
Basket Material	Polythene Coated Wire
Basket Size (mm) (WxDxH)	485mm x 180mm x 171mm
Interior Light	No
Fan	No

Spark Free	Yes
Defrost Setting	Auto (Off Cycle)
Controller Type	Digital, External
Bluetooth	Yes
Number Of Probes	1
Probe Type	Air
Ambient Temperature Range	10°C to 32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply (V)	230V
Insulation Thickness (mm)	45mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System Available	No
Optional Digital Lock	No
Other Door Option	No



For more information contact us on
 email: enquiries@medlinescientific.com
 Tel: +44 (0)1865 400321