

## Medline Analogue Hotplate (MS-100)

**Product Code:** MS100

**Brand:** Medline Scientific

Medline's Analogue Hotplate (MS-100) offers a large heating area of 180 x 180mm with safe controlled heating up to 400°C. Constructed from chemically resistant ceramic coated steel, this lab hotplate delivers excellent reliability with reproducible results at a sensible price. A standard 12.5mm screwed thread is included for the addition of a retort rod (not included).



## Specifications

<b>Description</b>	Hotplate
<b>Plate Area (mm)</b>	180 x 180
<b>Temperature Range (°C)</b>	Max. 400
<b>Dimensions (W x D x H)</b>	410 x 210 x 150mm
<b>Weight (kg)</b>	4

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)