

## Medline 100 to 1500rpm Analogue Stirrer (MS-200)

**Product Code:** MS200

**Brand:** Medline Scientific

Constructed from chemically resistant ceramic coated steel, Medline's 100 to 1500rpm Analogue Stirrer (MS-200) delivers excellent reliability with reproducible results at a sensible price. A standard 12.5mm screwed thread is included for the addition of a retort rod (not included).



## Specifications

<b>Description</b>	Stirrer
<b>Plate Area (mm)</b>	180 x 180
<b>Speed Range (rpm)</b>	100 to 1500
<b>Dimensions (W x D x H)</b>	410 x 210 x 150mm
<b>Weight (kg)</b>	4

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)