

Deltalab Sterile Amies Charcoal Swabs (Snappable Polystyrene & Viscose)

Product Code: 300285

Brand: Deltalab

Deltalab's Sterile Amies Charcoal Swabs (Snappable Polystyrene & Viscose) are a modification of Cary Blair media, which is at the same time a modification of Stuart media. It assures the viability of organisms such as *Trichomonas Vaginalis*, *Streptococcus Pyogenes*, and *Streptococcus Pneumoniae*.

Swabs are sterilised by radiation.



Specifications

Peel-Pack Dimensions (mm)	38 x 210
Material(s)	Snappable Polystyrene and Viscose (with charcoal)
Expiry Date	30 months after sterilisation date
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com