

## Medline 1.5ml Microcentrifuge Tubes

**Product Code:** 200400SP

**Brand:** Medline Scientific

Made from selected high-quality and biologically inert polypropylene, Medline's 1.5ml Microcentrifuge Tubes are suitable for use at 121°C/2.1 bar (cap not closed).

The leak-free caps on these graduated microcentrifuge tubes can be sealed efficiently and will stay closed during centrifugation with excellent visibility.

Come in packs of 1000.

### Features:

- Made from selected high-quality and biologically inert polypropylene.
- High resistance to temperatures and chemicals.
- Autoclavable at 121°C/2.1 bar (cap not closed).
- Place for markings on label area and flat cover.
- Perfect closure of leak free cap due to the design of the lid.
- Membrane can be punctured by sampling needles.
- Graduation.
- Round inside bottom for easy recovery of pellets or fluids.
- Easy closure (e.g. for use in chemical laboratories).



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)