

## Medline Printed Biohazard Specimen Bags

**Product Code:** MD66200

**Brand:** Medline Scientific

Made of high molecular and high-density polyethylene, Medline's Printed Biohazard Specimen Bags are translucent and can be used for the safe and convenient disposal of used pipettes, petri dishes, culture plates, membrane filters, and more.

### Specifications:

<b>Dimensions (w x l x h)</b>	140 x 150 x 200mm
<b>Thickness (mm)</b>	0.03
<b>Material</b>	Polyethylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)