

## Deltalab Sterilisation Indicator Tape (Ethylene Oxide)

**Product Code:** I91049

**Brand:** Deltalab

Deltalab's Sterilisation Indicator Tape (Ethylene Oxide) is self-adhesive and made of semi-creped paper. It also sticks on different surfaces such as plastic, metal, fabrics, paper, cardboard, etc. Will remain affixed for 18 months.

### Features:

- Overall thickness of 0.11 +/- 0.1 mm
- Indicator for sterilisation by ethylene oxide gas. The word GAS printed in brown changes to green once it has been sterilised as per the following specifications: exposition to 600 mg/l EO and 60%RH for minimum 20 minutes at 54 °C, and exposition to 600 mg/l EO and 60%RH for minimum 25 minutes at 37 °C

### Specifications:

<b>Type of Sterilisation</b>	Ethylene Oxide Gas
<b>Length (mm)</b>	50
<b>Width (mm)</b>	19
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)