

BECELEC 2

Electric Burner (3 Pin UK)

Product Code: BEC2UKSET



The **BECELEC 2 Electric Burner (3 PIN UK)** offer laboratories a unique safe replacement option to conventional gas Bunsen burners. The high power heater combined with a reflective internal stainless steel bowl produces a concentrated stream of heat comparable with that of existing gas-powered Bunsen burners.

Supplied with reflective internal stainless steel bowl, wire safety guard, retort stand, support rod, boss, and clamps. One litre gel flame optional.



Specifications:

Temperature Range (°C)	37 to 920
Diameter of Heating Area (mm)	Ø 130
Power	500 Watt
Safety Indicator	LED
Power Supply	220V AC
Plug Type	3 pin UK
Dimensions (w x d x h)	31 x 20 x 11cm
Weight (kg)	1.7
Accessories	Stainless steel plate, flame gel, support rod, boss head, and clamp

Phone: 0114 224 2257

Email: enquiriesUK@calibrescientific.com

Web: www.medlinescientific.co.uk