

# BECELEC 2

## Electric Burner (2 Pin EU)

Product Code: BEC2EXSET



The **BECELEC 2 Electric Burner (2 PIN EU)** offer laboratories a unique safe replacement option to conventional gas Bunsen burners. The high power heater combined with a reflective internal stainless steel bowl produces a concentrated stream of heat comparable with that of existing gas-powered Bunsen burners.

Supplied with reflective internal stainless steel bowl, wire safety guard, retort stand, support rod, boss, and clamps. One litre gel flame optional.



### Specifications:

Temperature Range (°C)	37 to 920
Diameter of Heating Area (mm)	Ø 130
Power	500 Watt
Safety Indicator	LED
Power Supply	220V AC
Plug Type	3 pin UK
Dimensions (w x d x h)	31 x 20 x 11cm
Weight (kg)	1.7
Accessories	Stainless steel plate, flame gel, support rod, boss head, and clamp

**Phone:** (+44) 0114 224 2257

**Email:** [enquiriesUK@calibrescientific.com](mailto:enquiriesUK@calibrescientific.com)

**Web:** [www.medlinescientific.co.uk](http://www.medlinescientific.co.uk)