

## Deltalab 120ml Aseptic Polypropylene Containers (57 x 73mm)

**Product Code:** 409702

**Brand:** Deltalab

Manufactured with ultra clear polypropylene, Deltalab's 120ml Aseptic Polypropylene Containers (57 x 73mm) feature a polyethylene leak proof red screw cap and a frosted panel on the side for easy writing.

Each sample container has moulded graduations in 20ml increments up to 100ml.

### Specifications

Recommended Volume (ml)	120
Description	Polypropylene
Sterility	Aseptic
Quantity	450

