

Deltalab Clear Containers with Screw Cap & No Label (200ml)

Product Code: 409902

Brand: Deltalab

Collect and store samples with Deltalab's Clear Containers with Screw Cap & No Label (200ml).

These clear sample containers have double moulded-in graduations (in 20ml increments up to 200ml and in 1oz increments up to 7oz) and feature a frosted panel on their side for easy writing.



Specifications:

Max. Capacity (ml)	240
Cap Colour	Red
Dimensions (mm)	62.2 x 115.5
Material(s)	Polypropylene (container) Polyethylene (cap)
Sterility	Not Sterile
Quantity	275

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com