

# Compound Microscopes – Lower Education

## Focus Monocular Compound Microscope

The Focus II monocular compound microscope is a school/college/university student microscope that comes with a built-in carrying handle. Economical with all the basic features required for educational use. Comfortable and reliable operation with modern design.



Comes with a **built-in carrying handle**

### Specifications

Cat. No.	1100.4440M	1100.4400M	1100.4000M
Optical Head	Monocular, inclined at 45° and 360° (rotatable)		
Siedentopf	Yes		
Dioptré	No		
Eyepieces	10x/18 mm widefield with pointer		
Turret Format	Triple		
Objectives	4x, 10x, 40x, 100x Achromatic	4x, 10x, 40x Achromatic	
Objective Type	Finite		
Condenser	Brightfield Abbe - NA 1.25 with iris diaphragm and filter holder. Rack and pinion adjustment for condenser		
Stage Size	120 x 120 mm with built-in condenser (N.A. 1.25)		
Focusing	Coaxial coarse and fine		
Illumination	20W halogen	0.5W LED	20W halogen
Package	One spare bulb, two spare fuses, dust cover, blue filter and instruction manual		