

Deltalab Blue Containers (60ml)

Product Code: 409511

Brand: Deltalab

Deltalab's Blue Containers (60ml) are ultra clear polypropylene, sterile containers with blue polyethylene leak proof screw caps.

Frosted panel on the side for easy writing.

Ideal for photosensitive samples or when it is required not to show the content.

Specifications:

Maximum Volume (ml)	60
Recommended Volume (ml)	50
Material	Polypropylene
Sterility	Sterile
Quantity	600



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com